




HPV

Sampling Solutions

HPV Sampling Solution-for Professionals

01 MSK-HPV001








Disposable Sampling Swab(s)

MS-VF4501 is a Flocked Swab used for specimen collection, and it has been designed for HPV sampling.

HPV Transport Medium

It is mainly composed of HBSS, Preservatives, Antibiotics, Ethylene Diamine Tetracetic Acid and Tritonx-100.

How to Use



HPV Self-sampling Solution

With or without Transport Mediums


- Perfect Design for Self-sampling
- Easy to Use and Improved Collection
- Save Time by not Going to the Hospital
- Solve Issue of Limited Collection Capacity



Support OEM / ODM

01

MSK-HPV002/HPV003



Disposable Sampling Swab(s)

MS-VF4503 is a Flocked Swab used for specimen collection, and it has been designed for HPV self-sampling.

HPV Transport Medium (2ml/5ml)

It is mainly composed of HBSS, Preservatives, Antibiotics, Ethylene Diamine Tetraacetic Acid and Tritonx-100.

02

MSK-HPV004



Simple Packing



Luxurious Packing

Disposable Sampling Swab(s)


MS-VF4503 is a Flocked Swab used for specimen collection, and it has been designed for HPV self-sampling.

HPV Transport Medium (3ml)

It is mainly composed of HBSS, Preservatives, Antibiotics, Ethylene Diamine Tetraacetic Acid and Tritonx-100.

03


MS-VF4504T



Self-sampling Swab packed by transport tube. After self-sampling, just put the swab into the transport tube.

04

MS-HSS001 / MSK-HSK001



Sampling Device



with Transport Medium

MS-HSS001 is a self-sampling device, and MSK-HSK001 is a self-sampling kit with HPV transport medium. It is not only easy to use, but also completely solves the issue of limited collection capacity while self-sampling.

Protect Your Privacy

Cervical Brushes and Kits

Cervical brushes are designed to collect cervical exfoliated cells, and they have been used for gynecological clinical diagnostics. The "tower" tip is mainly used for HPV screening, and the "escopa" tip is mainly used for TCT. MandeLab has developed the cervical brushes which exceed the industry standard according to the ergonomic design. The material of tip the is flexible and harmless, thus it would provide safe and comfortable sampling method. Cervical brushes collect more samples and ensure that sufficient samples can be released into the transport medium.



"Tower" Tip



"Escopa" Tip



Website: www.mandelab.com

E-mail: admin@mandelab.com

Tel: +86-755-28701995 +86-755-28701897

Add: 7F & 6F, Bld#5, 10th Industrial Zone, Tianliao Community,
Yutang Street, Guangming District, Shenzhen, 518132