

Environmental Sampling Kits

Environmental Sampling Kits have been used for collection, stabilization and transportation of microbial samples from a variety of surfaces. Environmental sampling kits contains 1 rayon swab and 1 tube with liquid Amies medium. The sampling system developed by MandeLab is highly suitable for automated & conventional processing. Customization is welcome.



Enviromental Sampling Kit Contains


Disposable Sampling Swab(s)

MS-VP006A is a Sampling Swab with Cotton tip & PP stick used for cell, virus and microbe specimen collection, and it also has been used for clinical diagnostics. Swabbing is one of the most widely used methods for microbiological examination of surface, so it has been used for bio-pharmaceutical industry, food-hygiene and cosmetics industry too. The head is made of medical grade cotton, it is high density and no easy to fall off, the PP handle is with suitable firmness and high resilience. The design is perfect and it is easy to use, the tip is harmless to micro-organisms, thus it would enhance the specimen collection.

Environmental Monitoring Transport Medium

The liquid transport medium is mainly composed of Potassium Dihydrogen Phosphate, Sodium Chloride, Peptone and Purified Water. It is used to preserve and transport the collected environmental samples, and only for in vitro analysis and testing, not for therapeutic use.

How to Use the Enviromental Sampling Kit?



Storage and Disposal

The optimum storage temperature is 2°C to 30 °C. Do not incubate, or freeze prior to use, improper storage will result in a loss of testing. Medical waste must be disposed of in compliance with the local legislation. Please take the appropriate precautions for infected material if necessary.

Product Pictures

